

Zea mays L. Embryo



GDS-80 No.: WGB090001
Barrel size: 4.5 mm
Pressure setting: 50 psi
Accessory: 3 cm target spacer
DNA concentration:
0.5 µg DNA/ 1.48 mg gold
Gold size: 1.6 µm
Reporter gene: GUS gene
GDS-80 setting: 4 rounds
Sample volume: 10 µL
Gas flow rate: 10~15 L/min

Sample description:

Fresh maize embryo, about 5 mm x 3 mm (L x W).



GDS-80 No.: WGB090001
Barrel size: 4.5 mm
Pressure setting: 50 psi
Accessory: 6 cm target spacer
DNA concentration:
0.5 µg DNA/ 1.48 mg gold
Gold size: 1.6 µm
Reporter gene: GUS gene
GDS-80 setting: 4 rounds
Sample volume: 10 µL
Gas flow rate: 10~15 L/min