

## Gower Ramsey Lip crest cells



GDS-80 No.: WGB10A-45008  
Barrel size: 4.5 mm  
Pressure setting: 40 psi  
Accessory: LC-10, 3 cm  
DNA concentration:  
1  $\mu$ g DNA/ 0.5 mg gold  
Gold size: 1.0  $\mu$ m  
Reporter gene: GUS gene  
GDS-80 setting: 5 rounds  
Sample volume: 10  $\mu$ L  
Gas flow rate: 10~15 L/min

## Sharry Baby Lip cells



GDS-80 No.: WGB10A-45008  
Barrel size: 4.5 mm  
Pressure setting: 40 psi  
Accessory: LC-10, 3 cm  
DNA concentration:  
1  $\mu$ g DNA/ 0.5 mg gold  
Gold size: 1.0  $\mu$ m  
Reporter gene: GUS gene  
GDS-80 setting: 5 rounds  
Sample volume: 10  $\mu$ L  
Gas flow rate: 10~15 L/min